



Lateral head immobilizer

USES

- | The lateral head immobilizer is used in combination with the spine board. Its main purpose is to prevent lateral movement of the head allowing to move the patient without the risk of aggravating cervical injuries.

DESIGN AND MATERIALS

- | This product has polyurethane foam stuffing that provides rigidity and volume to the object. This material can be removed through the side zippers of the coated polyester case, which facilitates cleaning and individual washing.
- | Fastening is done through polypropylene straps and acetal buckles, which allows the immobilizer to be secured to the spineboard.

