



Fastening straps

USES

This element was created to be used with spine boards, since they allow the immobilization of the patient for later transfer.

Its color coding ensures the correct use of the system. It allows for the regulation of the distance between anchorage points, as well as the adjustment on the body.

DESIGN AND MATERIALS

The straps are made of heavy 50 mm-wide polypropylene with Velcro to achieve a quick adjustment. It has a loop-type end, so it is easy to put it on and take it off.

These elements are translucent to X-rays, allowing the performance of diagnostic imaging without removing the immobilization elements, reducing the likelihood of posttraumatic injury.

