



Spider harness

USES

| A complimentary element to spinal boards that allows patient immobilization for safe transportation.

| The position of the straps and color codification ensures the correct use of the system. They allow to secure the patient in the areas of the shoulders, the chest, the abdomen, the hips and the legs, and provides the possibility of adjustment not only of the distance between the fastening points but also of the pressure applied to the body.

DESIGN AND MATERIALS

| It has 6 fastening straps differentiated by color in order to facilitate correct application and adjustment. They are made of heavy polypropylene with acetal buckles. These materials are translucent to X ray, allowing to perform diagnostic imaging without removing the immobilization elements, which reduces the likelihood of post-traumatic injury.

| This spider-type system can be used on any spinal board and on basket-type rescue stretchers.

ACCESSORIES

| Optional: a bag to store and transport it.

**Spider harness
placed on spinal board**

